



#### **General information**

**Description** The Y79 line was isolated by T.W. Reid and associates in January 1971 by explant culture of a primary tumor from

the right eye obtained immediately after enucleation. The donor had a strong maternal family history of retinoblastoma. Ultrastructural features including nuclear membrane infoldings, triple membrane structures, microtubules, large coated vesicles, centrioles, basal bodies and annulate lamellae were reportedly similar to

those of the original tumor.

Organism Human

**Tissue** Retina

**Disease** Retinoblastoma

**Synonyms** Y79, GM01232, GM01232E

#### **Characteristics**

Age 2,5 years

**Gender** Female

**Ethnicity** Caucasian

Morphology Multicellular clusters

Growth properties

**Clusters in Suspension** 

## **Identifiers / Biosafety / Citation**

**Citation** Y-79 (Cytion catalog number 300382)

Biosafety level

### **Expression / Mutation**

**Isoenzymes** PGM1, 1, G6PD, B, ES-D, 1, AK-1, 1, GLO-1, 2, Phenotype Frequency Product: 0.1373

Reverse transcriptase Negative



# Y-79 Cells | 300382

Karyotype	Hypertriploid, with abnormalities including dicentrics, breaks, pulverizations and minutes
Handling	
Culture Medium	RPMI 1640, w: 2.1 mM stable Glutamine, w: 2.0 g/L NaHCO3 (Cytion article number 820700a)
Medium supplements	Supplement the medium with 10% FBS
Subculturing	Maintain cultures by periodically adding or replacing the medium. Initiate cultures with a density of $2 \times 10^5$ cells/ml and keep the cell concentration within the range of $1 \times 10^5$ to $1 \times 10^6$ cells/ml for optimal growth.
Split ratio	A ratio of 1:2 to 1:4 is recommended
Seeding density	1x10^5 cells/ml
Fluid renewal	Every 3 days
Freeze medium	CM-1 (Cytion catalog number 800100) or CM-ACF (Cytion catalog number 806100)





#### Handling of cryopreserved cultures

- 1. Confirm that the vial remains deeply frozen upon delivery, as cells are shipped on dry ice to maintain optimal temperatures during transit.
- 2. Upon receipt, either store the cryovial immediately at temperatures below -150°C to ensure the preservation of cellular integrity, or proceed to step 3 if immediate culturing is required.
- 3. For immediate culturing, swiftly thaw the vial by immersing it in a 37°C water bath with clean water and an antimicrobial agent, agitating gently for 40-60 seconds until a small ice clump remains.
- 4. Perform all subsequent steps under sterile conditions in a flow hood, disinfecting the cryovial with 70% ethanol before opening.
- 5. Carefully open the disinfected vial and transfer the cell suspension into a 15 ml centrifuge tube containing 8 ml of room-temperature culture medium, mixing gently.
- 6. Centrifuge the mixture at 300 x g for 3 minutes to separate the cells and carefully discard the supernatant containing residual freezing medium. Optionally, skip centrifugation but remove any remaining freezing medium after 24 hours.
- 7. Gently resuspend the cell pellet in 10 ml of fresh culture medium. For adherent cells, divide the suspension between two T25 culture flasks; for suspension cultures, transfer all the medium into one T25 flask to promote effective cell interaction and growth.
- 8. Adhere to established subculture protocols for continued growth and maintenance of the cell line, ensuring reliable experimental outcomes.

# Quality control / Genetic profile / HLA

#### **Sterility**

Mycoplasma contamination is excluded using both PCR-based assays and luminescence-based mycoplasma detection methods.

To ensure there is no bacterial, fungal, or yeast contamination, cell cultures are subjected to daily visual inspections.



# Y-79 Cells | 300382

STR profile CSF1PO: 11,12

D13S317: 11,12
D16S539: 13,14
D5S818: 11,12
D7S820: 8,9
TH01: 6,9.3
TPOX: 8
vWA: 15,18
D3S1358: 15,16
D21S11: 30,32
D18S51: 13,16
Penta E: 13,18
Penta D: 12
D8S1179: 13,16
FGA: 22

**HLA alleles A\***: 02:01:01

**B\***: 01.01.1900 16:01, 02.01.1900 03:01

C\*: 03:04:01, 12:03:01

DRB1\*: 01:01:01, 13:02:01

DQA1\*: 01:01:01, 01:02:01

DQB1\*: 05:01:01, 06:04:01

DPB1\*: 03:01, 04:01

E: 01:01, 01:03

CLS Cell Lines Service GmbH | Dr.-Eckener-Str. 8 | 69214 Eppelheim | Germany Tel.: +49(0)6221 405780 | www.cytion.com | info@cytion.com