

Y-79 Cells | 300382

General information

Description The Y79 line was isolated by T.W. Reid and associates in January 1971 by explant culture of a primary tumor from the right eye obtained immediately after enucleation. The donor had a strong maternal family history of retinoblastoma. Ultrastructural features including nuclear membrane infoldings, triple membrane structures, microtubules, large coated vesicles, centrioles, basal bodies and annulate lamellae were reportedly similar to those of the original tumor.

Organism Human

Tissue Retina

Disease Retinoblastoma

Synonyms Y79, GM01232, GM01232E

Characteristics

Age 2,5 years

Gender Female

Ethnicity Caucasian

Morphology Multicellular clusters

Growth properties Clusters in Suspension

Regulatory Data

Citation Y-79 (Cytion catalog number 300382)

Biosafety level 1

NCBI_TaxID 9606

CellosaurusAccession CVCL_1893

Biomolecular Data

Y-79 Cells | 300382

Isoenzymes	PGM1, 1, G6PD, B, ES-D, 1, AK-1, 1, GLO-1, 2, Phenotype Frequency Product: 0.1373
Reverse transcriptase	Negative
Karyotype	Hypertriploid, with abnormalities including dicentrics, breaks, pulverizations and minutes

Handling

Culture Medium	RPMI 1640, w: 2.0 mM stable Glutamine, w: 2.0 g/L NaHCO3 (Cytion article number 820700a)
Subculturing	Maintain cultures by periodically adding or replacing the medium. Initiate cultures with a density of 2×10^5 cells/ml and keep the cell concentration within the range of 1×10^5 to 1×10^6 cells/ml for optimal growth.
Split ratio	A ratio of 1:2 to 1:4 is recommended
Seeding density	1×10^5 cells/ml
Fluid renewal	Every 3 days
Freeze medium	As a cryopreservation medium, use complete growth medium (including FBS) + 10% DMSO for adequate post-thaw viability, or CM-1 (Cytion catalog number 800100), which includes optimized osmoprotectants and metabolic stabilizers to enhance recovery and reduce cryo-induced stress.

Quality Control & Molecular Analysis

Sterility	<p>Mycoplasma contamination is excluded using both PCR-based assays and luminescence-based mycoplasma detection methods.</p> <p>To ensure there is no bacterial, fungal, or yeast contamination, cell cultures are subjected to daily visual inspections.</p>
------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Y-79 Cells | 300382

STR profile

CSF1PO: 11,12
D13S317: 11,12
D16S539: 13,14
D5S818: 11,12
D7S820: 8,9
TH01: 6,9.3
TPOX: 8
vWA: 15,18
D3S1358: 15,16
D21S11: 30,32
D18S51: 13,16
Penta E: 13,18
Penta D: 12
D8S1179: 13,16
FGA: 22

HLA alleles

A*: '02:01:01
B*: '40:01:02, '51:01:01
C*: '03:04:01, '12:03:01
DRB1*: '01:01:01, '13:02:01
DQA1*: '01:01:01, '01:02:01
DQB1*: '05:01:01, '06:04:01
DPB1*: '03:01:01, '04:01:01
E: '01:01, '01:03