



# Ba/F3 | 305224

**Culture Medium** RPMI 1640 2.0 2.0 / NaHCO<sub>3</sub> (820700a)

**Supplements** 5% FBS 10 IL-3

**Subculturing** 5, 5, 6 15 × 10

**Freeze medium** (FBS) + 10% DMSO

- Thawing and Culturing Cells**
1. ...
  2. ... -150
  3. ... 37
  4. ... 70%
  5. ... 15 ... 8
  6. ... 300 × 3
  7. ... 10
  8. ...

**Incubation Atmosphere** 37

**Flask Coating**

**Freezing Procedure**

Product sheet

Ba/F3 | 305224

**Shipping Conditions** [redacted] -78

**Storage Conditions** [redacted] -150 ° -196 [redacted]  
/ / HLA

**Sterility** [redacted] (PCR) [redacted]  
[redacted]