



3 | 305265

**Antigen expression** ICAM-1 + VCAM-1 + MAdCAM-1 + MAdCAM-1 +

**Viruses** Adenovirus (Ad2) (MPyV) T Adenovirus (PyMT)

**Culture Medium** DMEM 4.5 g/l 4 g/l 3.7 g/l NaHCO3 1.0 g/l (82)

**Supplements** 10% FBS


**Dissociation Reagent**

**Freeze medium** (10% FBS) + 10% DMSO

- Thawing and Culturing Cells**
1. ...
  2. ... -150 ...
  3. ... 37 ...
  4. ... 70% ...
  5. ... 15 ... 8 ...
  6. ... 300 × ... 3 ...
  7. ... 10 ...
  8. ...

Product sheet

3 | 305265

**Incubation Atmosphere** 37 


**Flask Coating** 

**Freezing Procedure** -78

**Shipping Conditions** -78

**Storage Conditions** -150 -196 

/ / HLA

**Sterility**  (PCR) 