

SCC-4 | 305384

Mutational profile TP53, p.Pro151Ser (c.451C>T)

XXXXXXXXXX

Culture Medium DMEM: DMEM:Ham's F12 (1:1) 3.1 / 2.5 15 (15)

Supplements 10 FBS 400 /

Dissociation Reagent

Subculturing PBS

Freeze medium (FBS) + 10% DMSO

Thawing and Culturing Cells

- 1. ...
- 2. ... -150
- 3. ... 37
- 4. ... 70%
- 5. ... 15 ... 8
- 6. ... 300 x 3
- 7. ... 10
- 8. ...

Incubation Atmosphere 37

Flask Coating

