

RS4:11 ~~XXXXXXXX~~ RS4:11 | 305360

~~XXXXXXXXXX~~ ~~XXXXXXXXXXXX~~ ~~XXXXXXXXXX~~

MSI-status ~~XX XXXXXXXX XX XXXXXXXX XXXXXXXX XXXX XXXXXXXX XXXXXXXX XXXXXXXX XXXXXXXX XXXX XXXXXXXX~~

~~XXXXXXXXXX~~

Culture Medium Alpha MEM ~~XX 2.0~~ ~~XXXXX XXXXXXXX XXXXXXXX XXXXXXXX XX XXXXXXXX XXXXXXXX XXXXXXXX XX 1.0~~ ~~XXXXX XXXXXXXX XXXXXXXX~~

Supplements ~~XXXX XXXXX XX 20% XX FBS XXXXXXXX XXXXXXXX XXXXX 20% XX XXXXXXXX XXXXXXXX XXXXXXXX XXXXXXXX~~

Split ratio ~~XXXXXX XXXXXXXX XXXXX XXXXXXXX XXXX 1:2 XX 1:4~~

Seeding density ~~XXXXXX XX ⁵XXXXXXXXXXXX-5 x 10~~

Fluid renewal 2 ~~XXXX 3 XXXXX XX XXXXXXXX~~

Freeze medium ~~XXXXXX XXXXXXXX XXXXXXXX XXXXXXXX XXXX XXXX XXXXX (XXXX XX XXX FBS) + 10% DMSO XX XXXX XXXXXXXX XXXX XXXXXXXX XXXXXXXX XXXX XXXXXXXX XXXXXXXX~~

- Thawing and Culturing Cells**
- ~~1. XXXXX XX XXXXX XXXXXXXX XXXXXXXX XXXXX XXXX XXX XXXXXXXX XXX XXX XXX XXXXXXXX XXX XXX XXX XXXXXXXX XXX XXXXX XXXXXXXX XXXXXXXX XXXXXXXX~~
 - ~~2. XXX XXXXXXXX XXX XXXXXXXX XXXXXXXX XXXXXXXX XXX XXXXXXX XX XXXXX XXXXXXX XXX -150 XXXXX XXXXXXX XXXXXXX XXXXXXX XXX XXXXXXX~~
 - ~~3. XXXXXXXX XXXXXXXX XX XXXXXXX XXXXXXXX XXXXXXX XX XXXXX XXXXXXX XX XXXXX XXXXX XXXXXXX XXXXXXX 37 XXXXX XXXXXXX XXXXX XXXXX XXXXXXX~~
 - ~~4. XXXXXXX XXXXX XXXXXXXX XXXXXXXX XX XXXXX XXXXXXX XX XXXXX XXXXXXX XX XXXXXXX XXXXXXXX XXXXXXXX XXXXXXX 70% XX XXXXXXXX XXXXXXX XXXXXXX~~
 - ~~5. XXXXX XXXXXXXX XXXXXXXX XXXXXXXX XXXXXXX XXXXX XXXXXXX XXX XXXXXXX XXX XXXXXXX XXX 15 XX XXXXXXX XXX 8 XX XX XXX XXXXXXXX XXXXXXX~~
 - ~~6. XXXXX XXXXXXXX XXXXXXXX XXXXXXXX XXX 300 x XX XXXXX 3 XXXXXXX XXXXX XXXXXXXX XXXXXXX XXXXXXX XX XXXXXXX XXXXXXXX XXXXXXX XXXXXXX XXXXXXX~~
 - ~~7. XXXXXXX XXXXXXX XXXXXXX XXXXXXX XXXXX XX 10 XX XX XXX XXXXXXX XXXXX XXXXXXXX XXXXXXXX XXXXXXXX XXX XXXXXXX XXXXXXX XXXXX XXXXXXX~~
 - ~~8. XXXXXXXX XXXXXXXX XXXXXXXX XXXXXXXX XXXXXXX XXX XXXXXXXX XXX XX XXXXXXX XXXXXXXX XXXXXXX XXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX~~

Incubation Atmosphere 37 ~~XXXXX XXXXXXXX XXX XX XXXXX~~

