

ND7/23 | 305520

[illegible]

ND7/23 | 305520

Seeding
density

1 - 3 × 10⁴ /

Freeze
medium

Freeze medium: 100% FBS (100% FBS) + 10% DMSO

Thawing and
Culturing Cells

1. Thaw cells quickly in a water bath at 37°C.
2. Add cells to a pre-warmed medium.
3. Incubate cells in a CO₂ incubator at 37°C.
4. Check cell morphology and confluency.
5. Pass cells when they reach 70-80% confluency.
6. Seed cells into new flasks.
7. Incubate cells in a CO₂ incubator at 37°C.
8. Check cell morphology and confluency.

Incubation
Atmosphere

37°C, 5% CO₂

Shipping
Conditions

Shipping conditions: -78°C

Storage
Conditions

Storage conditions: -150 to -196°C

Sterility

Sterility: PCR (PCR) and other methods.