

CHO-CXCR4 CHO-CXCR4 | 305411L

Product information

Description

CHO-CXCR4-Medium-High CHO (CHO) cells expressing human CXCR4. The cells are adapted to grow in DMEM/F12 medium supplemented with 10% fetal bovine serum (FBS) and 100 ng/ml insulin, transferrin, and selenium (ITS). The cells are maintained in DMEM/F12 medium supplemented with 10% FBS and 100 ng/ml ITS. The cells are adapted to grow in DMEM/F12 medium supplemented with 10% FBS and 100 ng/ml ITS. The cells are maintained in DMEM/F12 medium supplemented with 10% FBS and 100 ng/ml ITS.

Organism CHO

Tissue CHO

Synonyms CHO-CXCR4

Characteristics

Age 1-2 passages

Gender Male

Morphology Epithelial

Growth properties Adherent

Documentation

Citation CHO-CXCR4 (ATCC CRL-2919) (ATCC CRL-2919 305411MH)

Biosafety level 1

NCBI_TaxID 10029

GMO Status GMO-S1: This CHO line contains a recombinant construct enabling low-level expression of human CXCR4 for chemokine receptor studies. This classification applies only within Germany and may differ elsewhere.

Receptors expressed

Receptors expressed CXCR4 (CD184)

Product sheet

CHO-CXCR4 CHO-CXCR4 | 305411L

General Information

Culture Medium DMEM:DMEM:Ham's F12 (1:1) 3.1 µg/ml / 2.5 µg/ml 15 µg/ml InSCREENeX InSCREENeX INS-ME-1039

Supplements 5% FBS (G418-Sulfat) 0.5 µg/ml

Dissociation Reagent

Subculturing PBS

Fluid renewal 2-3

Post-Thaw Recovery 1:2 1:3 T25 70% 15 8

Freeze medium (FBS) + 10% DMSO

- Thawing and Culturing Cells**
1. Thaw cells in a 37°C water bath.
 2. Dilute cells into 10 ml of culture medium.
 3. Seed cells into a T25 flask.
 4. Allow cells to recover for 70% confluency.
 5. Harvest cells after 15-18 hours.
 6. Seed cells into a 300 cm² flask.
 7. Harvest cells after 10 days.
 8. Harvest cells after 10 days.

Incubation Atmosphere 37°C, 5% CO₂, humidified atmosphere.

Product sheet

CHO-CXCR4 CHO-CXCR4 | 305411L

Shipping Conditions

Cryopreserved cell lines are shipped on dry ice in validated, insulated packaging with sufficient refrigerant to maintain approximately -78 °C throughout transit. On receipt, inspect the container immediately and transfer vials without delay to appropriate storage.

Storage Conditions

For long-term preservation, place vials in vapor-phase liquid nitrogen at about -150 to -196 °C. Storage at -80 °C is acceptable only as a short interim step before transfer to liquid nitrogen.

HLA

Sterility

PCR