

XXXXXXXX SNU-423 SNU-423 | 305874

Antigen expression XXXXX XXXX B Rh ++

Mutational profile XXXX: TERT XXXXX c.1-124C>T (c.228C>T) (C228T) XXX XXXXX XXXXXXX=XX XXXXXXX. XXXXXXX TP53 XXXXX c.376-2A>G XXX XXXX

Karyotype XXXX XXXXXXX XXXXXXXXXX XXXXX XXXXXXXXXX = 79

XXXXXXXXXXXX

Culture Medium RPMI 1640 XX 2.0 XXXXX XXXXXXX XXXXXXXXXX XXXXXXX XX 2.0 XX/XXX NaHCO3 (XXXX XXXXXXX 820700a XX XXXXXXXXX)

Supplements XXXXX 10% XX FBS XXXXXXX XXXXXXXXXX XXX XXXXXXX

Dissociation Reagent XXXXXXX

Doubling time 72 XXXXX

Fluid renewal 2 XXXX 3 XXXXX XX XXXXXXXXXX

Freeze medium XXXXXXX XXXXXXX XXXXXXXXXX XXXXXXX XXX XXX XXXXX (XXXX XX XXX FBS) + 10% DMSO XX XXX XXXXXXX XXX XXXXXXX XXXXXXX XXX XXXXXXX XXXXXXX XXXXXXX

