

Product sheet

XXXXXXXX Lenti-X293T | 305820

Seeding density 2 $\times 10^4$ $\times \frac{1 \text{ ml}}{100 \mu\text{l}}$

Fluid renewal 2-3 times

Freeze medium XXXXX XXXXX XXXXXXXXXXXX XXXXXXXX XXX XXX XXXX + 10% DMSO XXXXXXXX XXX XXXXXXXX XXXXXXXX XXX XXXXXXXX XXXXXXXX

Thawing and Culturing Cells

1. XXXX XXX XXXXX XXXXXXXXXXXX XXXXXXXX XXXX XXX XXXXXXXXXXXX XXX XXX XXX XXXXXXXXXXXX XXX XXX XXX XXXXXXXXXXXX XXX XXXXXXXX XXX XXXXXXXX XXXXXXXXXXXX XXXXXXXX
2. XXX XXXXXXXXXXXX XXX XXXXXXXX XXXXXXXXXXXX XXXXXXXXXXXX XXX XXXXXXXX XXX XXXXXXXX XXX XXX -150 XXXXX XXXXXXXX XXXXXXXX XXXXXXXX XXX XXXXXXXX
3. XXXXXXXXXXXX XXXXXXXXXXXX XXX XXXXXXXX XXXXXXXXXXXX XXXXXXXX XXX XXXXXXXX XXX XXXXXXXX XXX XXXXXXXX XXXXXXXX 37 XXXXX XXXXXXXX XXXXXXXX XXXXXXXX
4. XXXXXXX XXXXX XXXXXXXXXXXX XXXXXXXX XXX XXXXXXX XXX XXXXXXX XXXXXXXXXXXX XXXXXXXXXXXX XXXXXXXXXXXX XXXXXXXX 70% XXX XXXXXXXXXXXX XXX XXXXXXXX
5. XXXXX XXXXXXXXXXXX XXXXXXXXXXXX XXXXXXXX XXXXXXX XXXX XXXXXXXXXXXX XXX XXXXXXX XXX XXXXXXX XXX 15 XXX XXXXXXX XXX 8 XXX XXX XXXXXXXXXXXX XXXXXXXX
6. XXX XXXXXXXX XXXXXXXXXXXX XXXXXXXX 200 \times XXX XXXXXXX 5 XXXXXXXXXXX XXXXXXX XXXXXXXXXXX XXX XXXXXXX XXXXXXXXXXX XXXXXXX XXXXXXX XXX XXXXXXXXXXXX
7. XXXXX XXXXXXXXXXX XXXXXXX XXX XXXXXXXXXXX XXX XXX XXXXXXXXXXX

Incubation Atmosphere 37 XXXXX XXXXXXXXXXXX XXX XXX XXXXXXX

Flask Coating XXX XXX

Shipping Conditions XXXX XXX XXXXXXXXXXX XXXXXXXXXXXX XXXXXXXXXXXX XXX XXX XXX XXX XXXXXXXXXXX XXXXXXXXXXX XXXXXXXXXXX XXX XXXXXXX XXXXXXXXXXX XXX XXXXXXX XXXXXXXXXXX -78

Storage Conditions XXXXXXX XXXXX XXXXXXXX XXX XXXXXXXXXXX XXX XXXXXXXXXXX XXXXXXX XXX XXXXXXX XXXXXXXXXXX XXX XXXXXXX XXXXXXX XXXXXXXXXXX XXX -150 \times -196 XXXXX XXXXXXX XXXXXXX

XXXXXXXX XXXXXXXX / XXXXXXXX XXXXXXXXXXXX / HLA