

# HEK293-CLDN6 | 305985

## HEK293-CLDN6

**Description**

HEK293-CLDN6 is a HEK293 cell line stably expressing CLDN6. HEK293-CLDN6 cells are derived from HEK293 cells (HEK293) and express CLDN6. HEK293-CLDN6 cells are used for studying CLDN6 function and for drug screening.

**Organism** Human

**Tissue** Kidney

## HEK293-CLDN6

**Age** 1-3

**Gender** Male

**Morphology** Epithelial

**Growth properties** Adherent

## HEK293-CLDN6

**Citation** HEK293-CLDN6 (HEK293-CLDN6) Cytion: 305985

**Biosafety level** 1

**NCBI\_TaxID** 9606

## HEK293-CLDN6

**Receptors expressed** CLDN6

## HEK293-CLDN6



Product sheet

HEK293-CLDN6 | 305985

**Storage Conditions**      Store at -150 °C to -196 °C

HLA

**Sterility**      Sterile (PCR)