

HEK293-CD20 | 305987

HEK293-CD20

Description

HEK293-CD20 is a HEK293 cell line expressing CD20. HEK293-CD20 cells are derived from HEK293 cells and express CD20. HEK293-CD20 cells are used for the production of recombinant antibodies. HEK293-CD20 cells are used for the production of recombinant antibodies. HEK293-CD20 cells are used for the production of recombinant antibodies.

Organism HEK293

Tissue HEK293

HEK293-CD20

Age HEK293

Gender HEK293

Morphology HEK293

Growth properties HEK293

HEK293-CD20

Citation HEK293-CD20 (HEK293-CD20) Cytion: 305987

Biosafety level 1

NCBI_TaxID 9606

HEK293-CD20

Receptors expressed HEK293

HEK293-CD20

Product sheet

HEK293-CD20 | 305987

Storage Conditions The cells should be stored at -150°C to -196°C in liquid nitrogen. The storage medium should be replaced every 6 months.

HEK293-CD20 / HEK293-CD20 / HLA

Sterility The cells are provided as a frozen cell suspension in a PCR tube. The cells are sterile and free of mycoplasma contamination. The cells are provided as a frozen cell suspension in a PCR tube. The cells are sterile and free of mycoplasma contamination.