

Product sheet

Vero E6 | 305008

Dissociation Reagent [REDACTED]

Doubling time 22 [REDACTED]

Subculturing [REDACTED] PBS [REDACTED]

Fluid renewal 2 [REDACTED] 3 [REDACTED]

Freeze medium [REDACTED] ([REDACTED] FBS) + 10% DMSO [REDACTED]

Thawing and Culturing Cells

1. [REDACTED]
2. [REDACTED] -150 [REDACTED]
3. [REDACTED] 37 [REDACTED]
4. [REDACTED] 70% [REDACTED]
5. [REDACTED] 15 [REDACTED] 8 [REDACTED]
6. [REDACTED] 300 × [REDACTED] 3 [REDACTED]
7. [REDACTED] 10 [REDACTED]
8. [REDACTED]

Incubation Atmosphere 37 [REDACTED]

Flask Coating [REDACTED]

Freezing Procedure [REDACTED] -78

Shipping Conditions [REDACTED] -78

Vero E6 | 305008

Storage Conditions -150 °C -196 °C

/ / HLA

Sterility (PCR)