

Product sheet

Yac-1 | 400383

General Information

Description	Yac-1 (Mo-MuLV) retrovirus
Organism	Yac-1
Tissue	Yac-1
Disease	Yac-1
Synonyms	YAC-1, YAC-1, YAC

Cell Line

Breed/Subspecies	Yac-1
Cell type	Yac-1
Growth properties	Yac-1

Identification

Citation	Yac-1 (400383)
Biosafety level	1
NCBI_TaxID	10090
CellSaurusAccession	CVCL_2244

Media and Supplements

Culture Medium

Culture Medium	RPMI 1640, 2.0, 2.0, NaHCO3 (820700a)
Supplements	10% FBS

Product sheet

Yac-1 | 400383

Subculturing 5 5 6 5 × 10

Fluid renewal 2 3

Freeze medium (FBS) + 10% DMSO

Thawing and Culturing Cells

- 1.
- 2. -150
- 3. 37
- 4. 70%
- 5. 15 8
- 6. 300 × 3
- 7. 10
- 8.

Incubation Atmosphere 37

Flask Coating

Freezing Procedure

Shipping Conditions

Storage Conditions -150 -196

HLA

Yac-1 | 400383

Sterility

Yac-1 (PCR) is a sterile, ready-to-use PCR master mix. It contains all the necessary components for PCR, including Taq polymerase, dNTPs, and buffer. The product is provided in a single-use format to ensure sterility and accuracy of results.