

TM4 | 305170

Antigen expression H-Y

Tumorigenic Yes

Culture Medium DMEM: DMEM:Ham's F12 (1:1) 3.1 / 2.5 15 (15)

Supplements 2.5% FBS 5%

Dissociation Reagent

Subculturing PBS


Fluid renewal 2 3

Freeze medium (FBS) + 10% DMSO


- Thawing and Culturing Cells**
- 1.
 2. -150
 3. 37
 4. 70%
 5. 15 8
 6. 300 × 3
 7. 10
 - 8.

TM4 | 305170

Incubation Atmosphere 37 

Flask Coating 

Freezing Procedure  -78

Shipping Conditions  -78

Storage Conditions  -150  -196 

/ / HLA

Sterility  (PCR) 
