

P3X63Ag8 | 305171

Description	Mouse myeloma cell line P3K27 (derived from P3K27 MOPC-21 fused with BALB/c).
Organism	Mouse
Tissue	Myeloma cell line derived from B cells
Synonyms	P3x63Ag8, P3-x63-Ag8, P3/x63-Ag8, P3/x63 Ag8, P3/x63/Ag8, P3-x63Ag8, P3x63 Ag8, P3x63 Ag8, P3 x 63Ag8, P3 x 63 Ag8, x63-Ag8, x63-AG8, x63-Ag8, P3x63, x63, GM03571
Breed/Subspecies	Swiss/129
Gender	Female
Morphology	Epithelial cells
Growth properties	Adherent
Citation	P3x63Ag8 (Cytion: 305171)
Biosafety level	1
NCBI_TaxID	10090
CellosaurusAccession	CVCL_3411

Protein expression	High
Antigen expression	2

P3X63Ag8 | 305171

[illegible]

P3X63Ag8 | 305171

Incubation Atmosphere 37 °C, 5% CO₂, humidified

Flask Coating Cell culture compatible

Shipping Conditions Dry ice, -78 °C

Storage Conditions -150 °C to -196 °C

Sterility Sterile, PCR compatible