

CHL | 305013

Dissociation Reagent [REDACTED]

Subculturing [REDACTED] PBS [REDACTED]

Fluid renewal 2 [REDACTED] 3 [REDACTED]

Freeze medium [REDACTED] ([REDACTED] FBS) + 10% DMSO [REDACTED]

- Thawing and Culturing Cells**
- [REDACTED]
 - [REDACTED] -150 [REDACTED]
 - [REDACTED] 37 [REDACTED]
 - [REDACTED] 70% [REDACTED]
 - [REDACTED] 15 [REDACTED] 8 [REDACTED]
 - [REDACTED] 300 × [REDACTED] 3 [REDACTED]
 - [REDACTED] 10 [REDACTED]
 - [REDACTED]

Incubation Atmosphere 37 [REDACTED]

Shipping Conditions [REDACTED] -78

Storage Conditions [REDACTED] -150 [REDACTED] -196 [REDACTED]

/ / HLA

