

Product sheet

HEK293T-KG-1a | 300234

Antigen expression	HLA A30A31B35Cw4
Isoenzymes	G6PD, B, PGM1, 1-2, PGM3, 0, ES-D, 1, Me-2, 1, AK-1, 0, GLO-1, 2
Viruses	EBVNA (EBNA):
Reverse transcriptase	
HEK293T-KG-1a	
Culture Medium	IMDM 4.5g/l, 4g/l, 25g/l (HEPES), 1.0g/l
Supplements	10% FBS
Doubling time	45h
Subculturing	5-6
Fluid renewal	3
Post-Thaw Recovery	24h
Freeze medium	10% DMSO + 10% FBS

