



**XXXXXXXX HROC300 T2 M1 HROC300 T2 M1 | 300866**

**XXXXXXXXXXXX**

**Culture Medium** DMEM: DMEM:Ham's F12 (1:1) 3.1 µg/ml / 2.5 µg/ml 15 µg/ml (15 µg/ml)

**Supplements** 10% FBS

**Dissociation Reagent** XXXXXXX

**Subculturing** XXXXXXX PBS XXXXXXX

**Fluid renewal** 3 5

**Freeze medium** XXXXXXX (XXX FBS) + 10% DMSO XXXXXXX

- Thawing and Culturing Cells**
1. XXXXXXX
  2. XXXXXXX -150 XXXXXXX
  3. XXXXXXX 37 XXXXXXX
  4. XXXXXXX 70% XXXXXXX
  5. XXXXXXX 15 XXXXXXX 8 XXXXXXX
  6. XXXXXXX 300 × XXXXXXX 3 XXXXXXX
  7. XXXXXXX 10 XXXXXXX
  8. XXXXXXX

**Incubation Atmosphere** 37 XXXXXXX

**Flask Coating** XXXXXXX

