

XXXXXXXX HROC300 T2 M1 HROC300 T2 M1 | 300866

XXXXXXXXXXXX

Culture Medium DMEM: DMEM:Ham's F12 (1:1) 3.1 µg/ml / 2.5 µg/ml 15 µg/ml (15 µg/ml)

Supplements 10% FBS

Dissociation Reagent XXXXXXX

Subculturing XXXXXXX PBS XXXXXXX

Fluid renewal 3 5

Freeze medium XXXXXXX (10% FBS) + 10% DMSO XXXXXXX

- Thawing and Culturing Cells**
1. XXXXXXX
 2. XXXXXXX -150 XXXXXXX
 3. XXXXXXX 37 XXXXXXX
 4. XXXXXXX 70% XXXXXXX
 5. XXXXXXX 15 XXXXXXX 8 XXXXXXX
 6. XXXXXXX 300 × 3 XXXXXXX
 7. XXXXXXX 10 XXXXXXX
 8. XXXXXXX

Incubation Atmosphere 37 XXXXXXX

Flask Coating XXXXXXX

