

Product sheet

Wnt-3A | 305184

Protein expression Wnt-3A

Wnt-3A

Culture Medium DMEM 4.5 g/l, 4 mM, 3.7 g/l, 1.0 mM (820)

Supplements 10% FBS, 0.4 mM G-418

Dissociation Reagent

Subculturing

Fluid renewal 2-3

Freeze medium (FBS) + 10% DMSO

- Thawing and Culturing Cells**
- 1.
 - 2.
 - 3.
 - 4.
 - 5.
 - 6.
 - 7.
 - 8.

Incubation Atmosphere 37

