

Product sheet

XXXXXXXX MH-S | 300487

Freezing Procedure XXXX XXXX XXXXXXXX XXXXXXXX XXXXXXXX XXXXXXXX XXX XXXX XXXX XXX XXXXXXX XXXXXXX XXXXXXXX XXX XXXXX XXXXX XXXXXXX XXX XXXXX XXXXXXX-78

Shipping Conditions XXXX XXXX XXXXXXX XXXXXXX XXXXXXXX XXXXXXXX XXX XXXX XXXX XXX XXXXXXX XXXXXXX XXXXXXXX XXX XXXXX XXXXX XXXXXXX XXX XXXXX XXXXXXX-78

Storage Conditions XXXXXXX XXXXX XXXXXXX XXX XXXXXXXX XXX XXXXXXXX XXXXX XXX XXXXXXX XXXXXXX XXX XXXXX XXXXXXX XXXXXXX XXX -150 ° -196 XXXXX XXXXXXX XXXXXXX

XXXXXXXX XXXXXXXX / XXXXXXX XXXXXXXX / HLA

Sterility XXXXXXX XXXXXXX XXXXXXXXXXXXXXXXXXXX XXXXXXXX XXX XXXXXXXXXXXXXXXXXXXX XXXXXXX XXXXXXXXXXXXXXXXXXXX XXXXXXXXXXXXXXXXXXXX (PCR) XXXXX XXXXXXX XXX
XXXXXXXX XXX XXXXX XXXXX XXXXXXX XXX XXXXX XXX XXXXXXX XXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX