

Product sheet

TE-1 | 305060

TE-1 | 305060

Description TE-1 is a cell line derived from a patient with a specific condition. It is characterized by its unique genetic profile and is used for research purposes. TE-1 cells are maintained in a specific culture medium and are known for their stability and reproducibility. TE-1 is a cell line derived from a patient with a specific condition. It is characterized by its unique genetic profile and is used for research purposes. TE-1 cells are maintained in a specific culture medium and are known for their stability and reproducibility.

Organism TE-1

Tissue TE-1

Disease TE-1

Synonyms TE1

TE-1 | 305060

Age 58

Gender TE-1

Ethnicity TE-1

Morphology TE-1

Growth properties TE-1

TE-1 | 305060

Citation TE-1 (TE-1 | 305060)

Biosafety level 1

NCBI_TaxID 9606

CellosaurusAccession CVCL_1759

TE-1 | 305060

TE-1 | 305060

TE-1

Culture Medium RPMI 1640 2.0 2.0 / NaHCO3 (820700a)

Supplements 10% FBS

Dissociation Reagent

Subculturing PBS

Fluid renewal 2-3

Freeze medium FBS + 10% DMSO

- Thawing and Culturing Cells
1. ...
2. ...
3. ...
4. ...
5. ...
6. ...
7. ...
8. ...

Incubation Atmosphere 37

Flask Coating

