

Product sheet

WEHI-164 | 400438

WEHI-164

Description WEHI-164 is a B cell lymphoma cell line derived from a BALB/c mouse. It is a non-adherent, non-tumorigenic cell line that grows in suspension. WEHI-164 cells are characterized by their high tumorigenicity and their ability to produce monoclonal antibodies. WEHI-164 cells are derived from a BALB/c mouse and are characterized by their high tumorigenicity and their ability to produce monoclonal antibodies.

Organism WEHI-164

Disease WEHI-164

Synonyms Wehi 164, Wehi164, Wehi 164 tc

WEHI-164

Breed/Subspecies WEHI-164

Morphology WEHI-164

Cell type WEHI-164

Growth properties WEHI-164

WEHI-164

Citation WEHI-164 (400438)

Biosafety level 1

NCBI_TaxID 10090

CellosaurusAccession CVCL_2251

WEHI-164

Tumorigenic WEHI-164 Balb/c

WEHI-164

