

Product sheet

WEHI-164 | 400438

WEHI-164

Description WEHI-164 is a murine B cell lymphoma cell line derived from a BALB/c mouse. It is a non-secreting, non-adherent cell line that grows in suspension. WEHI-164 cells are characterized by their ability to produce monoclonal antibodies. The cell line is maintained in culture with IL-4 and IL-6. WEHI-164 cells are used for the production of monoclonal antibodies against various antigens.

Organism Mouse

Disease B cell lymphoma

Synonyms Wehi 164, Wehi164, Wehi 164 tc

Characteristics

Breed/Subspecies BALB/c

Morphology Lymphoblastoid

Cell type B cell

Growth properties Suspension

Identification

Citation WEHI-164 (ATCC 400438)

Biosafety level 1

NCBI_TaxID 10090

CellosaurusAccession CVCL_2251

Characteristics

Tumorigenic Non-tumorigenic

Characteristics

