

Product sheet

TERA-C3H | 400360

General information

Description	Mouse, inbred strain, C3H, 4 weeks old
Organism	Mouse
Tissue	Whole body
Disease	None

Genetic background

Breed/Subspecies	C3H
Age	16 weeks
Gender	Male
Morphology	Normal
Growth properties	Normal

Identification and documentation

Citation	TERA-C3H (400360)
Biosafety level	1
NCBI_TaxID	10090
CellSaurusAccession	CVCL_AV91

Pathogen status and health

Tumorigenic	Not tumorigenic
Viruses	Not infected with MAP, M-Adenovirus, B.piliformis, etc.

XXXXXX TERA-C3H | 400360

XXXXXXXXXX

Culture Medium	DMEM 4.5 g/l, 4 g/l, 3.7 g/l, NaHCO3 1.0 g/l, (820 ml)
Supplements	10% FBS
Dissociation Reagent	
Subculturing	PBS
Fluid renewal	2
Freeze medium	(FBS) + 10% DMSO
Thawing and Culturing Cells	1. 37 °C 2. -150 °C 3. 37 °C 4. 70% 5. 15 6. 300 × 3 7. 10 8.
Incubation Atmosphere	37 °C
Flask Coating	

TERA-C3H | 400360

Freezing
Procedure

Das Produkt wird in einem Gefrierbehälter in einer Kühlkammer bei -78 °C gefroren.

Shipping
Conditions

Das Produkt wird in einem Gefrierbehälter in einer Kühlkammer bei -78 °C transportiert.

Storage
Conditions

Das Produkt wird in einem Gefrierbehälter in einer Kühlkammer bei -150 °C bis -196 °C gelagert.

Das Produkt ist in einem Gefrierbehälter verpackt.

Sterility

Das Produkt ist sterilisiert und kann für PCR-Analysen verwendet werden.
Das Produkt ist in einem Gefrierbehälter verpackt.