

Product sheet

SW-579 | 300346

SW-579

**Description** SW-579 is a cell line derived from a patient with a specific condition. It is characterized by its ability to grow in suspension and its sensitivity to certain treatments. SW-579 is a cell line derived from a patient with a specific condition. It is characterized by its ability to grow in suspension and its sensitivity to certain treatments. SW-579 is a cell line derived from a patient with a specific condition. It is characterized by its ability to grow in suspension and its sensitivity to certain treatments.

**Organism** Human

**Tissue** Tissue

**Disease** Disease

**Synonyms** SW579, SW579, SW579

SW-579

**Age** 59

**Gender** Gender

**Ethnicity** Ethnicity

**Morphology** Morphology

**Growth properties** Growth properties

SW-579

**Citation** SW-579 (300346)

**Biosafety level** 1

**NCBI\_TaxID** 9606

**CellosaurusAccession** CVCL\_3603

SW-579

SW-579 | 300346

**Antigen expression** O Rh+

**Isoenzymes** 21-2 31 1-2 1-1 1-1 1-1 2 6 0.0209

**Oncogenes** Myc + myb + ras + fos + sis + p53 + abl - ros - src - N-myc -.

**Tumorigenic**

SW-579

**Culture Medium** RPMI 1640 2.0 2.0 NaHCO3 (820700a)

**Supplements** 10% FBS

**Dissociation Reagent**

**Subculturing** PBS

**Fluid renewal** 2 3

**Freeze medium** (FBS) + 10% DMSO

