

Product sheet

RenCa-IL2 | 400322

Product Information

Description	RenCa-IL2 is a recombinant protein consisting of the extracellular domain of the interleukin-2 receptor (IL-2R) fused to the C-terminal of the RenCa protein. It is used for the treatment of RenCa in Balb/c mice. RenCa-IL2 is a recombinant protein consisting of the extracellular domain of the interleukin-2 receptor (IL-2R) fused to the C-terminal of the RenCa protein. It is used for the treatment of RenCa in Balb/c mice.
Organism	Mouse
Tissue	RenCa
Disease	RenCa
Synonyms	RENCA-IL-2

Animal Information

Breed/Subspecies	RenCa/IL2
Age	6 weeks
Gender	Male
Morphology	RenCa
Growth properties	RenCa

Identification

Citation	RenCa-IL2 (RenCa-IL2) 400322
Biosafety level	1
NCBI_TaxID	10090
CellosaurusAccession	CVCL_5944
GMO Status	GMO-S1: RenCa-IL2

Product sheet

RenCa-IL2 | 400322

XXXXXXXXXX XXXXXXXXXXXX XXXXXXXXXXXX

Tumorigenic XXXX XX XXXXXXXX XXXXXXXXXXXX

Products IL-2

XXXXXXXXXX

Culture Medium RPMI 1640 X 2.0 XXXXX XXXXX XXXXXXXXXXXX XXXXXXXX X 2.0 XX/XXX NaHCO3 (XXXX XXXXXXXX 820700a XX XXXXXXXX)

Supplements XX XXXXXXXX XXXXXXX X 10% XX XX FBS

Dissociation Reagent XXXXXXXX

Subculturing XX XXXXXXXX XXXXXXX XXXXXXX XX XXXXXXXX XXXXXXXXXXXX XXXXXXXX XXXXXXXXXXXX PBS XXXXX XXXXXXX XX XXXXXXXXXXXX XXXXXXXXXXXXXXXX XXXXXXXXXXXX XXXXXXX

Split ratio XXXXXXX XXXXXXXXXXXX XXXXX XXXXXXX XXX 1:4 X 1:8

Fluid renewal 2 XXX 3 XXXXX XX XXXXXXXX

Freeze medium XXXXXXX XXXXXXX XXXXXXXXXXXX XXXXXXX XXX XX XXXXX (XXXX XX XXX FBS) + 10% DMSO XX XXX XXXXXXX XX XXXXXXX XXXXXXX XXX XXXXXXX XXXXXXX

RenCa-IL2 | 400322

Thawing and Culturing Cells

1. [Redacted]
2. [Redacted]
3. [Redacted]
4. [Redacted]
5. [Redacted]
6. [Redacted]
7. [Redacted]
8. [Redacted]

Incubation Atmosphere 37 [Redacted]

Flask Coating [Redacted]

Freezing Procedure [Redacted]

Shipping Conditions [Redacted]

Storage Conditions [Redacted]

[Redacted] / [Redacted] / HLA

Sterility [Redacted]

[Redacted] STRAmelogenin: x/y