



**XXXXXXXX RCC-MH | 300237**

**Receptors expressed** HLA-A2 XXXXX

**Protein expression** IL8

**Mutational profile** IL8 RS1126647 3-UTR Wt

**XXXXXXXXXXXX**

**Culture Medium** RPMI 1640 X 2.0 XXXXX XXXXXX XXXXXXXXXXXX XXXXXXXX X 2.0 X/X NaHCO3 (XXXX XXXXXXXX 820700a XX XXXXXXXX)

**Supplements** XX XXXXXXX XXXXXX X 10X XX X FBS

**Dissociation Reagent** XXXXXXX

**Subculturing** XX XXXXXXX XXXXXX XXXXXXX XX XXXXXXX XXXXXXXXXXX XXXXXXX XXXXXXXXXXXX PBS XXXX XXXXXX XX XXXXXXXXXXXX XXXXXXXXXXXX XXXXXXXXXXXX XXXXXXXXXXXX

**Split ratio** XXXX XXXXXX 1:2 XXX 1:3

**Fluid renewal** 1 XXX 2 XXXX XX XXXXXXXX

**Freeze medium** XXXXXXX XXXXXXX XXXXXXXXXXXX XXXXXXX XXX XXX XXXX (XXXX XX XXX FBS) + 10% DMSO XX XXX XXXXXXX XXX XXXXXXX XXXXXXX XXXXXXXXXXXX XXXXXXXXXXXX XXXXXXXXXXXX



**XXXXXXXX RCC-MH | 300237**

---

**XXXXXXXX XXXXXXXXXXXX STRAmelogenin: x**

**CSF1PO:** 11,12

**D13S317:** 11

**D16S539:** 12

**D5S818:** 9,1

**D7S820:** 8,12

**TH01:** 7

**TPOX:** 8

**vWA:** 15,18

**D3S1358:** 17

**D21S11:** 29,3

**D18S51:** 16

**Penta E:** 5,12

**Penta D:** 12,13

**D8S1179:** 9:15

**FGA:** 22

**PEZ6:** HROG17