





CCF-STTG1 | 300388

Thawing and Culturing Cells

1. Thaw the cells in a water bath at 37°C. Transfer the cells to a centrifuge tube and centrifuge at 300 x g for 3 minutes. Remove the supernatant and resuspend the cells in 10 ml of complete medium. Seed the cells into a T75 flask and incubate at 37°C in 5% CO2. The next day, replace the medium with fresh complete medium.
2. Once the cells have reached confluence, trypsinize and seed into a 96-well plate at 70% confluence. Incubate at 37°C in 5% CO2.
3. Harvest the cells at 15 days post-infection. Collect the supernatant and store at -150°C. The cells should be harvested at 8 days post-infection.
4. Harvest the cells at 10 days post-infection. Collect the supernatant and store at -150°C. The cells should be harvested at 8 days post-infection.
5. Harvest the cells at 10 days post-infection. Collect the supernatant and store at -150°C. The cells should be harvested at 8 days post-infection.
6. Harvest the cells at 10 days post-infection. Collect the supernatant and store at -150°C. The cells should be harvested at 8 days post-infection.
7. Harvest the cells at 10 days post-infection. Collect the supernatant and store at -150°C. The cells should be harvested at 8 days post-infection.
8. Harvest the cells at 10 days post-infection. Collect the supernatant and store at -150°C. The cells should be harvested at 8 days post-infection.

Incubation Atmosphere

37°C, 5% CO2

Flask Coating

Flask coating is not required for this cell line.

Freezing Procedure

Freeze the cells in a freezing medium and store at -80°C.

Shipping Conditions

Ship the cells at -80°C.

Storage Conditions

Store the cells at -150°C to -196°C.

CCF-STTG1 / HLA

Sterility

The cells are sterile and free of mycoplasmas. PCR testing is performed to ensure sterility.

XXXXXXXX CCF-STTG1 | 300388

---

XXXXXXXX HLA

**A\***: '01:01:01

**B\***: '08:01:01, '37:01:01

**C\***: '06:02:01, '07:01:01

**DRB1\***: '07:01:01, '13:02:01

**DQA1\***: '01:02:01, '02:01:01

**DQB1\***: '03:03:02, '06:04:01

**DPB1\***: '04:01:01

**E**: '01:01:01