

SKW-3 | 300343

Antigen expression

CD2+ CD3- CD3- CD4+ CD8⁺ Thy-1

Products

2 (2 2)

Culture Medium

RPMI 1640 2.0 2.0 NaHCO₃ (820700a)

Supplements

10% FBS

Doubling time

30

Subculturing

5 5 6 5 x 10

Post-Thaw Recovery

1 x 10⁵

Freeze medium

10% DMSO + 10% FBS

Thawing and Culturing Cells

1. 10
2. -150
3. 37
4. 70%
5. 15 8
6. 300 x 3
7. 10
8. 10

████████ SKW-3 | 300343

████████ HLA

A*: '11:01:01, '30:01:01
B*: '35:01:01, '44:02:01
C*: '04:01:01, '05:01:01
DRB1*: '01:03:01, '04:01:01
DQA1*: '01:01:01, '03:03:01
DQB1*: '03:01, '05:01
DPB1*: '04:01:01, '04:02:01
E: '01:01:01