

5 ( ) : 3.0 / :  
: 2.0 : 2.2 /  
| 820200a

McCoy's 5A Medium 3.0 2.2 820200a

McCoy's 5A (3.0) 3.0 2.2 NaHCO

- 7.2 +/- 0.02 20-25
- 37
- 6 8

|                   |                      |       |
|-------------------|----------------------|-------|
|                   | /                    |       |
| McCoy's 5A Medium | 3.0 2.2 820200a      | 132   |
|                   | 37                   | 97    |
|                   | 6 8                  | 400   |
|                   | 6,460                |       |
|                   | 504                  |       |
| McCoy's 5A Medium | D(+)-3.0 2.2 820200a | 3,000 |
|                   | 0                    |       |
|                   | 600                  |       |
|                   | 11                   |       |

Product sheet

5 ( ) : 3.0 / :  
: 2.0 : 2.2 /  
| 820200a

|                      |                                       |       |
|----------------------|---------------------------------------|-------|
| XXXXXXXXXXXXXXXXXXXX | L-XXXXXXXXXXXX                        | 13    |
|                      | L-XXXXXXXXXXXX × HCl                  | 42.14 |
|                      | L-XXXXXXXXXXXX × H <sub>2</sub>       | 45    |
|                      | XXX L-XXXXXXXXXXXX                    | 19.9  |
|                      | L-XXXXXXXXXXXX × $\frac{HCl}{2}$ × H  | 31.7  |
|                      | L-XXXXXXXXXXXX (XXXXXX)               | 219.1 |
|                      | XXX L-XXXXXXXXXXXX                    | 22    |
|                      | XXXXXXXXXXXX                          | 7     |
|                      | L-XXXXXXXXXXXX × OHCl × H             | 20.9  |
|                      | L-XXXXXXXXXXXXXXXXXXXX                | 19.6  |
|                      | L-XXXXXXXXXXXX                        | 39    |
|                      | L-XXXXXX                              | 39.36 |
|                      | L-XXXXXX × HCl                        | 36.54 |
|                      | L-XXXXXXXXXXXX                        | 14.9  |
|                      | L-XXXXXX XXXXXXXX                     | 16.5  |
|                      | L-XXXXXXXXXXXX                        | 17.27 |
|                      | L-XXXXXX                              | 26.28 |
|                      | L-XXXXXXXXXXXX                        | 17.87 |
|                      | L-XXXXXXXXXXXX                        | 3     |
|                      | XXX XXXXXXX XXXXXXXXXX L-XXXXXXXXXXXX | 26.10 |

Product sheet

5 ( ) : 3.0 / :  
: 2.0 : 2.2 /  
| 820200a

|               |                                |      |
|---------------|--------------------------------|------|
|               | L- <span></span>               | 17.5 |
| <span></span> | <span></span> p- <span></span> | 1    |
|               | <span></span> <span></span>    | 0    |
|               | D(+)- <span></span>            | 0    |
|               | D- <span></span> <span></span> | 0    |
|               | <span></span> <span></span>    | 5    |
|               | <span></span> <span></span>    | 10   |
|               | <span></span> - <span></span>  | 36   |
|               | <span></span>                  | 0    |
|               | <span></span> <span></span>    | 0    |
|               | <span></span> <span></span>    | 0    |
|               | <span></span> <span></span>    | 0.5  |
|               | <span></span>                  | 0    |
|               | <span></span> <span></span>    | 0    |
|               | <span></span> <span></span>    | 2    |