

HEP2 BJAB | 302006

HEP2 **HEP2** **HEP2**

Antigen expression +CD10+ +CD19+ +CD20+ +CD21(+) +CD22+ +CD23- +CD24- +CD32+ +CD37+ +CD38+ +CD39- +CD40+ +CD54+ +CD72+ +CD73- +CD75+ +CD77+ +CD81+ +CD82+ +CD83+ +CD84+ +CD86

Karyotype 46,XX,XX,XX,XX,XX

HEP2

Culture Medium RPMI 1640, 2.0 mM L-glutamine, 2.0 mM NaHCO₃ (820700a)

Supplements 20% FBS, 10 mM beta-mercaptoethanol

Subculturing 1:5 to 1:6

Seeding density 3 x 10⁵ cells

Fluid renewal 3-5 days

Post-Thaw Recovery 48 hours

Freeze medium 10% FBS + 10% DMSO

XXXXXXXXBJAB | 302006

XXXXXXXX HLA

A*: 01:01:83, '02:01:01
B*: '13:02:01, '35:01:01
C*: '04:01:01, '06:02:01
DRB1*: '12:01:01, '13:02:01
DQA1*: '01:02:01, '05:05:01
DQB1*: '03:01, '06:04:01
DPB1*: '04:02:01G
E: '01:01, '01:03