



# GIMEN | 300179

**Culture Medium** DMEM 4.5 g/l, 4 g/l, 3.7 g/l, NaHCO<sub>3</sub> 1.0 g/l (82)

**Supplements** 10% FBS

**Dissociation Reagent**

**Doubling time** 25

**Subculturing** PBS

**Seeding density** 2 × 3 × 4


**Fluid renewal** 2 × 3

**Freeze medium** (FBS) + 10% DMSO

- Thawing and Culturing Cells**
1. ...
  2. ... -150 ...
  3. ... 37 ...
  4. ... 70% ...
  5. ... 15 ... 8 ...
  6. ... 300 × 3 ...
  7. ... 10 ...
  8. ...

GIMEN | 300179

**Incubation Atmosphere** 37 

**Flask Coating** 



**Freezing Procedure** -78

**Shipping Conditions** -78

**Storage Conditions** -150 -196 

/ / HLA

**Sterility**  (PCR)   


- HLA**
- A\*:** '02:01:01, '30:01:01
- B\*:** '13:02:01, '18:01:01
- C\*:** '06:02:01, '07:01:09
- DRB1\*:** '04:03:01, '07:01:01
- DQA1\*:** '02:01:01, '03:01:01
- DQB1\*:** '02:02:01, '03:02:01
- DPB1\*:** '02:01:02
- E:** '01:01:01'01:01:01'01:xx