

GIMEN | 300179

Culture Medium DMEM 4.5 g/l, 4 g/l, 3.7 g/l, NaHCO₃ 1.0 g/l (82)

Supplements 10% FBS

Dissociation Reagent

Doubling time 25

Subculturing PBS

Seeding density 2 × 10⁴ - 3 × 10⁵

Fluid renewal 2 × 3

Freeze medium (10% FBS) + 10% DMSO

- Thawing and Culturing Cells**
1. ...
 2. ... -150 ...
 3. ... 37 ...
 4. ... 70% ...
 5. ... 15 ... 8 ...
 6. ... 300 × ... 3 ...
 7. ... 10 ...
 8. ...

GIMEN | 300179

Incubation Atmosphere

37

Flask Coating

Freezing Procedure

-78

Shipping Conditions

-78

Storage Conditions

-150 -196

/ / HLA

Sterility

(PCR)

HLA

- A*: '02:01:01, '30:01:01
- B*: '13:02:01, '18:01:01
- C*: '06:02:01, '07:01:09
- DRB1*: '04:03:01, '07:01:01
- DQA1*: '02:01:01, '03:01:01
- DQB1*: '02:02:01, '03:02:01
- DPB1*: '02:01:02
- E: '01:01:01 '01:01:01 '01:xx