

XXXXXXXX XXXXXXXX | 300366

Thawing and Culturing Cells

1. XXXXX XX XXXXX XXXXXXXXXXXX XXXXXXX XXXX XXX XXXXXXXXXXXX XXX XXX XXX XXXXXXXXXXXX XXX XXX XXX XXXXXXXXXXXX XXX XXXXXXX XXXXXXXXXXXX XXXXXXXXXXXX XXXXXXX
2. XXX XXXXXXXXXXXX XXX XXXXXXX XXXXXXXXXXXX XXXXXXXXXXXX XXX XXXXXXX XX XXXXX XXXXXXX XXX -150 XXXXX XXXXXXX XXXXXXX XXXXXXX XXX XXXXXXX
3. XXXXXXXXXXXX XXXXXXXXXXXX XX XXXXXXXXXXXX XXXXXXXXXXXX XXXXXXX XX XXXXX XXXXXXX XX XXXXXXX XXXXXXX XXXXXXX XXXXXXX 37 XXXXXXX XXXXXXX XXXXX XXXXXXX XXX
4. XXXXXXX XXXXX XXXXXXXXXXXX XXXXXXXXXXXX XX XXXXX XXXXXXX XX XXXXX XXXXXXX XX XXXXXXX XXXXXXXXXXXX XXXXXXXXXXXX XXXXXXX 70% XX XXXXXXXXXXXX XXXXXXX
5. XXXXX XXXXXXXXXXXX XXXXXXXXXXXX XXXXXXXXXXXX XXXXXXX XXXXX XXXXXXXXXXXX XXX XXXXXXX XXX XXXXXXX XXX 15 XX XXXXXXX XXX 8 XX XX XXX XXXXXXXXXXXX XXX
6. XXXXX XXXXXXXXXXXX XXXXXXXXXXXX XXX 300 × XX XXXXX 3 XXXXXXX XXXXX XXXXXXXXXXXX XXXXXXX XXXXXXXXXXXX XX XXXXXXXXXXXX XXXXXXXXXXXX XXXXX XXXXXXX XXXXXXX
7. XXXXXXX XXXXXXX XXXXXXX XXXXXXXXXXX XXXXX XX 10 XX XX XXX XXXXXXX XXXXX. XXXXXXXXXXXX XXXXXXXXXXXX XXXXXXXXXXXX XX XXXXXXXXXXXX XXXXXXXXXXXX XXX XXXXXXX
8. XXXXXXXXXXXX XXXXXXXXXXXX XXXXXXXXXXXX XXXXXXXXXXXX XXX XXXXXXXXXXXX XXX XX XXXXXXXXXXXX XXXXXXXXXXXX XXXXXXX XXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX

Incubation Atmosphere 37 XXXXX XXXXXXXXXXXX XXX XXX XXXXX
XXXXX

Flask Coating XX XXXX

Freezing Procedure XXXX XXX XXXXXXXXXXXX XXXXXXXXXXXX XXXXXXXXXXXX XXX XXXXX XXX XX XXXXXXX XXXXXXXXXXXX XXXXXXXXXXXX XX XXXXXXX XXXXX XXXXXXXXXXXX XXX XXXXXXX XXXXXXX -78

Shipping Conditions XXXX XXX XXXXXXXXXXXX XXXXXXXXXXXX XXXXXXXXXXXX XXX XXXXX XXX XX XXXXXXX XXXXXXXXXXXX XXXXXXXXXXXX XX XXXXXXX XXXXX XXXXXXXXXXXX XXX XXXXXXX XXXXXXX -78

Storage Conditions XXXXXXX XXXXX XXXXXXXXXXXX XX XXXXXXXXXXXX XXX XXXXXXXXXXXX XXXXX XX XXXXXXX XXXXXXXXXXXX XXX XXXXXXX XXXXXXX XXXXXXXXXXXX XXX -150 X -196 XXXXXXX XXXXXXX XXXXXXX

XXXXXXXX XXXXXXXX / XXXXXXX XXXXXXXXXXXX / HLA

Sterility XXXXXXXXXXX XXXXXXX XXXXXXXXXXXXXXXXXXXX XXXXXXXXXXXX XX XX XXXXXXXXXXXX XXXXXXXXXXXX XXX XXXXXXX XXXXXXXXXXXXXXXXXXXX XXXXXXXXXXXX (PCR) XXXXX XXXXXXX XXX
XXXXXXXX XXX XXXXX XXXXX XXXXXXXXXXX XX XXXXX XX XXXXXXX XXX XXXXXXX XXXXXXX XXXXXXXXXXX XXXXXXX XXXXXXXXXXX XXXXXXXXXXX

XXXXXXXX XXXXXXXX | 300366

XXXXXXXX HLA

A*: '03:01:01, '74:01:01

B*: '53:01:01, '58:01:01

C*: '04:01:01

DRB1*: '11:02:01, '15:03:01

DQA1*: '01:02:01, '05:05:01

DQB1*: '03:19:01, '06:02:01

DPB1*: '01:01:01, '02:01:02

E: '01:01, '01:03