

Product sheet

XXXXXXXX MKN-74 | 300490

Storage Conditions XXXXX XXXX XXXXXX XX XXXXXXXXXXXX XX XXXXXXXXXXXX XXXX XX XXXXXX XXXXXXXXXXX XXX XXXX XXXXXX XXXXXXXXXXX XXXX -150 X -196 XXXXX XXXXXX XXXXXX

XXXXXXXX XXXXXXXXXXX / XXXXXXX XXXXXXXXXXX / HLA

Sterility XXXXXXXXXXX XXXXXXX XXXXXXXXXXXXXXXXXXXX XXXXXXXXXXXX XX XX XXXXXXXXXXXX XXXXXXXXXXX XXX XXXXXX XXXXXXXXXXXXXXXXXXXX XXXXXXXXXXXX (PCR) XXXXX XXXXXX XX
XXXXXXXX XX XXXX XXXX XXXXXX XX XXXXX XX XXXXXX XXX XXXXXX XXXXXX XXXXXXXXXXX XXXXXX XXXXXXXXXXX XXXXXXXXXXX