

P3X63Ag8.653 P3X63Ag8.653 | 400118

XXXXXXXXXX

Culture Medium RPMI 1640 2.0 2.0 NaHCO3 (820700a)

Supplements 10% FBS

Dissociation Reagent

Subculturing 15 min PBS

Seeding density 2 x 10

Fluid renewal

Post-Thaw Recovery 48

Freeze medium 10% DMSO

- 1. ...
- 2. ... -150
- 3. ... 37
- 4. ... 70%
- 5. ... 15 ... 8
- 6. ... 300 x 3
- 7. ... 10
- 8. ...

