

BRL | 305193

Dissociation Reagent [REDACTED]

Subculturing [REDACTED] PBS [REDACTED]

Fluid renewal 2 [REDACTED] 3 [REDACTED]

Freeze medium [REDACTED] ([REDACTED] FBS) + 10% DMSO [REDACTED]

- Thawing and Culturing Cells**
1. [REDACTED]
 2. [REDACTED] -150 [REDACTED]
 3. [REDACTED] 37 [REDACTED]
 4. [REDACTED] 70% [REDACTED]
 5. [REDACTED] 15 [REDACTED] 8 [REDACTED]
 6. [REDACTED] 300 × [REDACTED] 3 [REDACTED]
 7. [REDACTED] 10 [REDACTED]
 8. [REDACTED]

Incubation Atmosphere 37 [REDACTED]

Flask Coating [REDACTED]

Freezing Procedure [REDACTED] -78

Shipping Conditions [REDACTED] -78

BRL | 305193

**Storage
Conditions**

Store at -150 °C to -196 °C. Do not freeze.

/ / HLA

Sterility

Not sterile. For PCR applications, the product must be treated with a DNAse treatment kit (PCR) before use.

Do not use if the product is contaminated or if the packaging is damaged.