

Product sheet

XXXXXXXXXX L-138 | 400384

Storage Conditions XXXXX XXXX XXXXXXXX XX XXXXXXXXXXXX XX XXXXXXXXXXXX XXXX XX XXXXXXX XXXXXXXX XXX XXXXX XXXXXXX XXXXXXXX XXXX -150 X -196 XXXXX XXXXXXX XXXXXXX

XXXXXXXXX XXXXXXXX / XXXXXXXX XXXXXXXX / HLA

Sterility XXXXXXXX XXXXXXXX XXXXXXXXXXXXXXXXXXXX XXXXXXXXXXXX XX XX XXXXXXXXXXXX XXXXXXXX XXX XXXXXXX XXXXXXXXXXXXXXXXXXXX XXXXXXXXXXXX (PCR) XXXXX XXXXXXX XX
XXXXXXXXX XXX XXXX XXXX XXXXXXX XX XXXXX XX XXXXXXX XXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX