

XXXXXXXXL-WRN | 300641

XXXXXXXXXX XXXXXXXXXXXX XXXXXXXXXXXX

Protein expression Wnt-3A, R-spondin, noggin

XXXXXXXXXX

Culture Medium DMEM 4.5 g/l / 4 g/l NaHCO_3 1.0 g/l (82)

Supplements 10% FBS

Dissociation Reagent

Subculturing

Freeze medium

Thawing and Culturing Cells

1. XXXXX
2. XXXXX -150
3. XXXXX 37
4. XXXXX 70%
5. XXXXX 15 8
6. XXXXX 300 x 3
7. XXXXX 10
8. XXXXX

Incubation Atmosphere 37

