

Product sheet

U-343 MG | 300365

U-343 MG

Description U-343 MG is a recombinant protein that is secreted by the cells. It is a dimeric protein with a molecular weight of approximately 100 kDa. The protein is highly stable and can be stored at -20°C for several months. It is used for the study of the function of the protein and for the development of diagnostic and therapeutic applications.

Organism *Escherichia coli*

Tissue *Escherichia coli*

Disease *Escherichia coli*

Synonyms U-343mg, u-343-mg, u-343-mg, u-343mg, u-343, u-343, u343, 343 mg, 343mg

U-343 MG

Age 54 days

Gender *Escherichia coli*

Ethnicity *Escherichia coli*

Morphology *Escherichia coli*

Growth properties *Escherichia coli*

U-343 MG

Citation U-343 MG (U-343 MG | 300365)

Biosafety level 1

NCBI_TaxID 9606

CellosaurusAccession CVCL_S471

U-343 MG

XXXXXXXXU-343 MG | 300365

XXXXXXXX HLA

A*: '02:01:01, '03:01:01

B*: '07:02:01, '47:01:01

C*: '06:02:01, '07:02:01

DRB1*: '04:05:01, '15:01:01

DQA1*: '01:02:01, '03:03:01

DQB1*: '03:01, '06:02

DPB1*: '04:01:01

E: '01:01:01