

HSC-T6 Cells | 305199

General information

Description

The HSC-T6 cell line is a well-characterized hepatic stellate cell line derived from adult rat liver tissue. These cells play a critical role in liver physiology and pathology, particularly in the processes of liver fibrosis and cirrhosis. Hepatic stellate cells are responsible for the storage of vitamin A in lipid droplets under normal physiological conditions. Upon liver injury, they transdifferentiate into myofibroblast-like cells, which secrete extracellular matrix proteins, contributing to the fibrotic response. The HSC-T6 cell line has been extensively utilized as a model to study these mechanisms due to its ability to mimic the in vivo behavior of activated hepatic stellate cells.

HSC-T6 cells express key markers such as α -smooth muscle actin (α -SMA), glial fibrillary acidic protein (GFAP), and desmin, which are indicative of their myofibroblastic phenotype. These cells also exhibit significant proliferative capacity and are responsive to various cytokines and growth factors, making them an invaluable tool for investigating the signaling pathways involved in hepatic fibrosis. Researchers have employed HSC-T6 cells to explore therapeutic targets and interventions aimed at mitigating fibrosis and promoting liver regeneration. The availability of this cell line has thus facilitated significant advancements in the understanding of liver disease and the development of potential treatments.

Organism Rat

Tissue Liver

Synonyms HSCT6

Characteristics

Breed/Subspecies Sprague Dawley

Age Adult

Gender Male

Morphology Epithelial

Growth properties Adherent

Regulatory Data

Citation HSC-T6 (Cytion catalog number 305199)

Biosafety level 1

HSC-T6 Cells | 305199

NCBI_TaxID	10116
------------	-------

CellosaurusAccession	CVCL_0315
----------------------	-----------

Biomolecular Data

Handling

Culture Medium	DMEM, w: 4.5 g/L Glucose, w: 4 mM L-Glutamine, w: 3.7 g/L NaHCO ₃ , w: 1.0 mM Sodium pyruvate (Cytion article number 820300a)
-----------------------	--

Subculturing	Remove the old medium from the adherent cells and wash them with PBS that lacks calcium and magnesium. For T25 flasks, use 3-5 ml of PBS, and for T75 flasks, use 5-10 ml. Then, cover the cells completely with Accutase, using 1-2 ml for T25 flasks and 2.5 ml for T75 flasks. Let the cells incubate at room temperature for 8-10 minutes to detach them. After incubation, gently mix the cells with 10 ml of medium to resuspend them, then centrifuge at 300xg for 3 minutes. Discard the supernatant, resuspend the cells in fresh medium, and transfer them into new flasks that already contain fresh medium.
---------------------	---

Split ratio	1:2 to 1:4
--------------------	------------

Fluid renewal	2 to 3 times per week
----------------------	-----------------------

Freeze medium	As a cryopreservation medium, use complete growth medium (including FBS) + 10% DMSO for adequate post-thaw viability, or CM-1 (Cytion catalog number 800100), which includes optimized osmoprotectants and metabolic stabilizers to enhance recovery and reduce cryo-induced stress.
----------------------	--

Quality Control & Molecular Analysis

Sterility	<p>Mycoplasma contamination is excluded using both PCR-based assays and luminescence-based mycoplasma detection methods.</p> <p>To ensure there is no bacterial, fungal, or yeast contamination, cell cultures are subjected to daily visual inspections.</p>
------------------	---