

NCI-H2452 Cells | 300391

General information

Description

The NCI-H2452 cell line is a human malignant pleural mesothelioma cell line, which was derived from the pleura of a patient with mesothelioma. It is frequently used in research focused on understanding mesothelioma pathophysiology and developing new therapeutic approaches. Like other mesothelioma cell lines, NCI-H2452 is associated with exposure to asbestos fibers, a well-established risk factor for mesothelioma. Studies involving NCI-H2452 have highlighted its utility in exploring mechanisms of disease progression and response to various therapies, particularly gene therapies and viral oncolysis approaches.

NCI-H2452 cells express Coxsackie and adenovirus receptor (CAR) and CD46, which makes them suitable candidates for adenovirus-based gene therapy studies. In research investigating oncolytic virotherapy, both adenovirus type 5 (Ad5) and a fiber-modified variant (Ad5F35) have been tested on NCI-H2452 cells. These adenoviruses replicate selectively within tumor cells, inducing oncolysis in a viral particle-dependent manner. It was found that both Ad5 and Ad5F35 displayed similar efficacy in inducing cell death in NCI-H2452 cells, supporting their potential in gene therapy for malignant mesothelioma.

In addition to its role in oncolytic virotherapy, NCI-H2452 cells have been used to study tumor angiogenesis, a key factor in mesothelioma progression. NCI-H2452 expresses progranulin (PGRN) and granulin-like proteins, which have been identified as novel angiogenic factors that operate independently of the VEGF pathway. This VEGF-independent angiogenesis is crucial, as it offers alternative therapeutic targets in cases where anti-VEGF therapies like bevacizumab fail to improve patient outcomes. Research indicates that these granulins contribute significantly to the formation of new blood vessels, which supports tumor growth and may be involved in the resistance to certain treatments.

Organism	Human
Tissue	Lung
Disease	Pleural biphasic mesothelioma
Synonyms	NCI-H2452, H-2452, NCIH2452

Characteristics

Age	Adult
Gender	Male
Ethnicity	European
Morphology	Epithelial
Growth properties	Adherent

NCI-H2452 Cells | 300391

Regulatory Data

Citation	NCI-H2452 (Cytion catalog number 300391)
-----------------	--

NCBI_TaxID	9606
-------------------	------

CellosaurusAccession	CVCL_1553
-----------------------------	-----------

Biomolecular Data

Handling

Culture Medium	RPMI 1640, w: 2.0 mM stable Glutamine, w: 2.0 g/L NaHCO ₃ (Cytion article number 820700a)
-----------------------	--

Supplements	Supplement the medium with 10% FBS
--------------------	------------------------------------

Dissociation Reagent	Accutase
-----------------------------	----------

Subculturing	Remove the old medium from the adherent cells and wash them with PBS that lacks calcium and magnesium. For T25 flasks, use 3-5 ml of PBS, and for T75 flasks, use 5-10 ml. Then, cover the cells completely with Accutase, using 1-2 ml for T25 flasks and 2.5 ml for T75 flasks. Let the cells incubate at room temperature for 8-10 minutes to detach them. After incubation, gently mix the cells with 10 ml of medium to resuspend them, then centrifuge at 300xg for 3 minutes. Discard the supernatant, resuspend the cells in fresh medium, and transfer them into new flasks that already contain fresh medium.
---------------------	---

Freeze medium	As a cryopreservation medium, we use complete growth medium (including FBS) + 10% DMSO for adequate post-thaw viability, or CM-1 (Cytion catalog number 800100), which includes optimized osmoprotectants and metabolic stabilizers to enhance recovery and reduce cryo-induced stress.
----------------------	---

NCI-H2452 Cells | 300391

Thawing and Culturing Cells

1. Confirm that the vial remains deeply frozen upon delivery, as cells are shipped on dry ice to maintain optimal temperatures during transit.
2. Upon receipt, either store the cryovial immediately at temperatures below -150°C to ensure the preservation of cellular integrity, or proceed to step 3 if immediate culturing is required.
3. For immediate culturing, swiftly thaw the vial by immersing it in a 37°C water bath with clean water and an antimicrobial agent, agitating gently for 40-60 seconds until a small ice clump remains.
4. Perform all subsequent steps under sterile conditions in a flow hood, disinfecting the cryovial with 70% ethanol before opening.
5. Carefully open the disinfected vial and transfer the cell suspension into a 15 ml centrifuge tube containing 8 ml of room-temperature culture medium, mixing gently.
6. Centrifuge the mixture at $300 \times g$ for 3 minutes to separate the cells and carefully discard the supernatant containing residual freezing medium.
7. Gently resuspend the cell pellet in 10 ml of fresh culture medium. For adherent cells, divide the suspension between two T25 culture flasks; for suspension cultures, transfer all the medium into one T25 flask to promote effective cell interaction and growth.
8. Adhere to established subculture protocols for continued growth and maintenance of the cell line, ensuring reliable experimental outcomes.

Incubation Atmosphere

37°C , 5% CO_2 , humidified atmosphere.

Shipping Conditions

Cryopreserved cell lines are shipped on dry ice in validated, insulated packaging with sufficient refrigerant to maintain approximately -78°C throughout transit. On receipt, inspect the container immediately and transfer vials without delay to appropriate storage.

Storage Conditions

For long-term preservation, place vials in vapor-phase liquid nitrogen at about -150 to -196°C . Storage at -80°C is acceptable only as a short interim step before transfer to liquid nitrogen.

Quality Control & Molecular Analysis

NCI-H2452 Cells | 300391

Sterility

Mycoplasma contamination is excluded using both PCR-based assays and luminescence-based mycoplasma detection methods.

To ensure there is no bacterial, fungal, or yeast contamination, cell cultures are subjected to daily visual inspections.