

**HEK293T L6 | 305231**

**HEK293T L6**

|             |                                                                                                                                                                                                                                                                                                                                                                             |
|-------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description | HEK293T L6 cells are a derivative of HEK293 cells, which are a derivative of HEK293 cells. HEK293 cells are a derivative of HEK293 cells, which are a derivative of HEK293 cells. HEK293T L6 cells are a derivative of HEK293T cells, which are a derivative of HEK293T cells. HEK293T L6 cells are a derivative of HEK293T cells, which are a derivative of HEK293T cells. |
| Organism    | HEK293T                                                                                                                                                                                                                                                                                                                                                                     |
| Tissue      | HEK293T                                                                                                                                                                                                                                                                                                                                                                     |
| Synonyms    | L-6, L-6 HEK293T                                                                                                                                                                                                                                                                                                                                                            |

**HEK293T L6**

|                   |         |
|-------------------|---------|
| Age               | HEK293T |
| Gender            | HEK293T |
| Cell type         | HEK293T |
| Growth properties | HEK293T |

**HEK293T L6**

|                      |                            |
|----------------------|----------------------------|
| Citation             | L6 (HEK293T) Cytion 305231 |
| Biosafety level      | 1                          |
| NCBI_TaxID           | 10116                      |
| CellosaurusAccession | CVCL_0385                  |

**HEK293T L6**

|                    |         |
|--------------------|---------|
| Protein expression | HEK293T |
|--------------------|---------|

**HEK293T L6**

[illegible]

**HLA L6 | 305231**

**Shipping  
Conditions**

HLA L6 is shipped in a dry ice container. The container must be kept at -78°C.

**Storage  
Conditions**

HLA L6 should be stored at -150 to -196 °C. The storage time is unlimited.

**HLA L6 / HLA L6 / HLA**

**Sterility**

HLA L6 is sterile and ready for use. PCR products are not required. HLA L6 is not a virus and does not contain any infectious agents. HLA L6 is not a virus and does not contain any infectious agents.