

HEK293T SNU-423 | 305874

Antigen expression HEK293T B; Rh +

Mutational profile HEK293T: TERT, HEK293T, c.1-124C>T (c.228C>T) (C228T), HEK293T, HEK293T=HEK293T HEK293T, TP53, HEK293T, c.376-2A>G, HEK293T, HEK293T

Karyotype HEK293T; HEK293T = 79

HEK293T

Culture Medium RPMI 1640, w: 2.0 mM HEK293T HEK293T, w: 2.0 g/L NaHCO3 (HEK293T HEK293T HEK293T Cytion 820700a)

Supplements HEK293T HEK293T 10% FBS HEK293T HEK293T

Dissociation Reagent HEK293T

Doubling time 72 HEK293T

Fluid renewal 2 HEK293T 3 HEK293T HEK293T

Freeze medium HEK293T HEK293T HEK293T, HEK293T HEK293T HEK293T HEK293T (HEK293T FBS) + 10% DMSO HEK293T HEK293T HEK293T HEK293T HEK293T HEK293T, HEK293T

