

Product sheet

XXXXXXXX D-XXXXXXXX (XXXXXXXXXXXXXXXX), 250 XXXX/XXXX | 831105

XXXXXXXX XXXXXXXX

XXXX 50 °F XXXXXXX +2°C XXX+8°C XXXXXXX: XXXXX XXXXXXX

XXXX D-(XXXXXXXXXXXXXXXX) XXXXX XX XXXXX XXXXX D-(+)-XXXXX, XXXXX XXXXX XXXXXXXXXXX XXXX. XXXXXXX XXXXX XXXXXXXXXXX XXXXXXXXXXX XXXXXXX

XXXXXXXXXXXXXXXXXXXXXXXXXXXX

XXXXX XXXXXXX XX XXXXXXX XXXXXXX XXXX XXXXXXX XXXXXXX XXXX XXXXXXX XX XX XXXXXXX XXXXXXXXXXX XXXX XXXXXXXXXXX. XXXXX XXXXX XXXXX XXXXXXX, XXXXXXX XXXX

XXXXX XXXXXXX

XXXXXXXX XXXXXXX XXXXXXX XXXX 250 XXXX/XXXX XXXXXXX, XX XXXXXXX XXXXXXX XXXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX. XXX XXXXXXX XXXXXXXXXXX XXX XXXXXXX XXXXXXX XX