

HEPES (1M) | 831106

100 g, +2°C to +8°C, sealed, sterile

HEPES (1M) is N-2-hydroxyethylpiperazine-N-2-ethanesulfonic acid, a zwitterionic buffer.

HEPES is a good buffer for most biological applications, especially for pH 7.0 to 8.0. It is a good buffer for CO₂ incubation. 10–25 mM HEPES is recommended for most applications. It is a good buffer for cell culture media.

1 M HEPES is a good buffer for most biological applications, especially for pH 7.0 to 8.0. It is a good buffer for CO₂ incubation. 10–25 mM HEPES is recommended for most applications. It is a good buffer for cell culture media.