





Product sheet

AAV-293 | 305127

Freezing Procedure

AAV-293 is a lentiviral vector. It is produced in HEK293 cells. The virus is purified by ultracentrifugation and stored at -78°C.

Shipping Conditions

AAV-293 is shipped in a dry ice container. The container is labeled with biohazard symbols and the text "Dry Ice -78°C".

Storage Conditions

AAV-293 is stored at -150 to -196°C in a liquid nitrogen storage container. The container is labeled with biohazard symbols and the text "Liquid Nitrogen -150 to -196°C".

AAV-293 / AAV-293 / HLA

Sterility

AAV-293 is produced in HEK293 cells. The virus is purified by ultracentrifugation and stored at -78°C. The virus is sterile and free of endotoxins. The virus is tested for sterility by PCR and endotoxin testing.