



















Product sheet

5637  | 300105

Freeze medium  (FBS) + 10% DMSO 

Thawing and Culturing Cells

- 
-  -150°C 
-  37 
-  70% 
-  15  8 
-  -300 x g  3 
-  -10  T2 
- 

Incubation Atmosphere 37°C, 5% CO₂ 

Flask Coating 

Freezing Procedure  -78°C 



Shipping Conditions  -78°C 

Storage Conditions  -150  196 

 /  / HLA

5637  | 300105

Sterility

 HLA

A*: '11:01:01, '68:02:01
B*: 15:03:01, 55:02:01
C*: '01:02:01, '02:10:01
DRB1*: '01:02:01, '09:01:02G
DQA1*: '01:01:02, '03:02:01
DQB1*: '03:03:02, '05:01:01
DPB1*: '05:01:01G, '13:01:01G
E: 01:03:02